



Drum Tipper


The drum tipper is used to empty drums with pharmaceutical tablets in a fixed hopper. The drum tipper lifts the drum, forms a dust-tight seal between the drum and the equipment, tips the drum and discharges bulk tablets. The standard unit is intended for dumping pharmaceutical tablets and pellets from drums up to 45 kg but can be redesigned according to your specifications..

TECHNICO

PHARMACEUTICAL ENGINEERING B.V.

Kerkstraat 34
4265JE Genderen
The Netherlands

Telephone: +31(0)416-352571
Mobile: +31(0)6-20049480
E-mail: cas.dernison@technicopharma.com



If the tablets in the hopper are to be kept in a protective atmospheric condition, the hopper needs to be provided with special valves in the lid to fill it. The Drum Tipper will control those valves. A closed hopper with pneumatic valves can also be supplied by Technico Pharma.

The Technico Pharma Drum Tipper can be ordered with one or two tipping units, suitable for different drum brands like Curtec, Lobo, Remcon.



Check the website or scan the QR code for pictures of many machine variants and a demonstration movie!

<https://www.technicopharma.com/machines/drum-tipper-dumper/4>



WWW.TECHNICOPHARMA.COM